HumaPette

Gebrauchsanleitung · Operating Manual
Mode d'emploi · Instrucciones de manejo

Human

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Table of Contents

Scope of supply	000
1 Terms of use	29
1.1 Hazard levels	29
1.2 Format	29
2 Safety regulations	30
2.1 General safety regulations	30
2.2 Purpose	31
2.3 Limitations of use	31
2.4 Operating limitations	31
2.5 Operating exclusions	31
3 Functions and controls	32
4 Pipetting	33
5 Checking the volume	36
6 Accuracy table	39
7 Calibration	40
8 Cleaning and disinfection	42
8.1 Autoclaving	42
8.2 UV sterilization	42
8.3 PE filter	42
9 Maintenance	43
9.1 Disassembling/cleaning (up to 1,000 µl)	43
9.2 Disassembling/cleaning (5 ml/10 ml)	45
10 Troubleshooting	46
11 Product markings	46
12 Ordering Information	47
12.1 Order info/accessories	47
12.2 Spare parts	48
12.2.1 Up to 1,000 µl	48
12.2.2 5 ml/10 ml	49
12.3 Additional accessories	50

13 Repairs	50
13.1 Sending for repair	50
14 Warranty	50
15 Disposal	51

1 Terms of use

- Please carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.

1.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

1.2 Format

Format	Meaning	Icon	Meaning
1. Task	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The HumaPette can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and observe this operating manual before using the device.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves. When working with infectious or hazardous samples, the standard laboratory rules and precautions must be adhered to.
3. Follow the instructions given by the reagent manufacturer.
4. Use the device only for pipetting liquids within the defined limitations and restrictions of use. Comply with the operating exclusions; see *Operating exclusions, p. 31*. In case of doubt, contact the manufacturer or dealer.
5. Always perform work in a manner that does not endanger users or other people. Avoid splattering. Use only suitable vessels.
6. Avoid touching the tip opening when working with aggressive media.
7. Never use force.
8. Use only original spare parts. Do not make any technical modifications. Do not disassemble the device further than described in the user manual!
9. Always check that the device is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks), stop pipetting immediately and refer to the section "Troubleshooting"; see *Troubleshooting, p. 46*. Contact the manufacturer, if necessary.

2.2 Purpose

This is an air displacement pipette for pipetting aqueous solutions of medium density and low to medium viscosity.

2.3 Limitations of use

This instrument is intended for pipetting samples, within the following limitations:

- Operating temperature of instrument and reagent should be between +15 °C and +40 °C (59 °F to 104 °F) (other temperatures upon request)
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s

2.4 Operating limitations


Viscous and wetting liquids may compromise volumetric accuracy. Volumetric accuracy may also be affected when pipetting liquids whose temperature deviates from the ambient temperature by more than ± 1 °C/ ± 1.8 °F.

2.5 Operating exclusions

The user is responsible for checking the compatibility of the device with the intended application. The device cannot be used:

- for liquids that corrode polypropylene
- for liquids that corrode polycarbonate (viewing window)
- for liquids with very high steam pressure
- for liquids that corrode FKM and polyether ether ketone (PEEK)

3 Functions and controls

1 Pipetting button

2 Volume-setting wheel

3 Cover

4 Finger rest

5 Tip cone

6 Pipetting shaft

7 Hand grip

8 Volume display

9 Tip ejection button

Volume display

The numbers in the display are read from top to bottom; the white dash corresponds to the decimal point.

Tip ejection button

The serial number is found behind the tip ejection button.

4 Pipetting

1. Inserting a tip

NOTICE


- 5 ml and 10 ml devices should only be used with a built-in PE filter; see *PE filter*, p. 42.
- Perfect analysis results can only be achieved by using quality tips. We recommend pipette tips from HUMAN. For additional information, refer to the Accuracy table.
- Pipette tips are disposable products!

a.

- a. Use the correct tips, in accordance with the volume range or color code! Make sure that the tips are firmly in place and leak tight.

2. Setting the volume


- a. Turn the volume-setting wheel to select the desired volume. In doing so, turn the adjustment wheel steadily, avoiding abrupt turning motions.


3. Aspirating a sample

NOTICE

The ISO 8655 standard requires that pipette tips are pre-wetted once before the actual pipetting procedure.

- a. Press the pipetting button until first resistance is felt.
- b. Hold the device vertically and immerse the tip in the liquid.


- c. Allow the pipetting button to steadily move back to its original position. Leave the tip immersed in the liquid for a few seconds, so that the set volume is aspirated completely. This is especially important when pipetting viscous media and when using pipettes with large volumes.

Volume range	Immersion depth [mm]	Wait time [s]
0.1 μl - 1 μl	1 - 2	1
> 1 μl - 100 μl	2 - 3	1
> 100 μl - 1000 μl	2 - 4	1
> 1,000 μl	3 - 6	3


NOTICE

Do not lay the instrument down when the tip is filled; this can cause the medium to flow into the instrument and contaminate it! The instrument should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

4. Dispensing a sample

- a. Place the pipette tip against the vessel wall. Hold the pipette at an angle of 30-45° to the vessel wall.
- b. Press the pipetting button at a uniform speed until the first resistance is felt and hold it. To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density media.


- c. Completely empty the tip by over-stroking: press the pipetting button until the second resistance is felt.
- d. While doing this, wipe the pipette tip against the vessel wall.
- e. Remove the pipette tip from the vessel wall and allow the pipetting button to move back to its original position.

5. Ejecting a tip

NOTICE

The device should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

- a. Hold the pipette shaft over a suitable disposal bin and press the tip ejection button all the way down.

5 Checking the volume

We recommend testing the device every 3 to 12 months, depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655, Part 6.

1. Setting the nominal volume

Set the maximum specified device volume (refer to page 8 for procedure).

2. Conditioning the pipette

Condition the pipette before testing by aspirating and dispensing the test liquid (H₂O distilled) with a pipette tip five times.

3. Performing the test

NOTICE

In accordance with DIN EN ISO 8655-2, a tip change is recommended after each individual measurement. An exception to this rule can be made, according to DAkkS guideline DKD-R8-1.

- a. Aspirate the test liquid and pipette into the weighing vessel.
- b. Weigh the pipetted amount with an analysis scale. (Please refer to the user manual of the scale manufacturer.)
- c. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- d. At least 10 pipetting series and weighings in 3 volume ranges (100 %, 50 %, 10 %) are recommended.

Calculation (for nominal volume) x_i = Weighing results n = Number of weighings Z = Correction factor(e.g. 1.0029 $\mu\text{l}/\text{mg}$ at 20 °C, 1013 hPa)

$$\text{Mean} \quad \bar{x} = \frac{\sum x_i}{n}$$

$$\text{Mean volume} \quad \bar{V} = \bar{x} \cdot z$$

$$\text{Accuracy}^* \quad A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

 V_0 = Nominal volume

$$\text{CV}\% = \frac{\quad}{\bar{V}}$$

$$\text{Standard deviation}^* \quad s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$


*) = Calculation for accuracy (A%) coefficient of variation (CV%): A% and CV% are calculated using the formulas of statistical quality control.

6 Accuracy table

Adjustable model

Volume range [µl]	Partial volume [µl]	A* ≤ ± %	CV* ≤ %	Sub steps [µl]	Recommended tip type [µl]
0.5 - 10	10	1	0.5	0.01	0.5 - 20
	5	1.6	1		
	1	7	4		
2 - 20	20	0.8	0.4	0.02	2 - 200
	10	1.2	0.7		
	2	5	2		
10 - 100	100	0.6	0.2	0.1	2 - 200
	50	0.8	0.4		
	10	3	1		
20 - 200	200	0.6	0.2	0.2	2 - 200
	100	0.8	0.3		
	20	3	0.6		
100 - 1000	1000	0.6	0.2	1	50 - 1000
	500	0.8	0.3		
	100	3	0.6		
500 - 5000	5000	0.6	0.2	5	500 - 5000
	2500	0.8	0.3		
	500	3	0.6		
1000 - 10,000	10000	0.6	0.2	10	1000 - 10000
	5000	0.8	0.3		
	1000	3	0.6		

*A = Accuracy, CV = Coefficient of Variation

Final test values based on the nominal volume (= max. volume) printed on the instrument and the specified partial volumes at the same temperature (20 °C/68 °F) of the instrument, surroundings and distilled water, in accordance with DIN EN ISO 8655.

Fixed-volume model

Volume range [μl]	A* ≤ ± %	CV* ≤ %	Recommended tip type [μl]
10	1	0.5	0.5 - 20
20	0.8	0.4	2 - 200
25	0.8	0.4	2 - 200
50	0.8	0.4	2 - 200
100	0.6	0.2	2 - 200
200	0.6	0.2	2 - 200
500	0.6	0.2	50 - 1000
1000	0.6	0.2	50 - 1000


*A = Accuracy, CV = Coefficient of Variation

7 Calibration

The device is permanently calibrated for aqueous solutions. The pipette can be calibrated if it is operating inaccurately or must be adjusted to work with solutions of varying density and viscosity or with specially-shaped pipette tips.


- a. Perform a volume check and determine the actual value; see *Checking the volume*, p. 36.

b.

b. Removing the cover: Push the clamp forward, gently lift and pull backward.

d.


d. Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

e.

e. Set the adjustment value: Make corrections by turning the volume-setting wheel in +/- direction. A volume check is recommended after every adjustment.

f.

f. Slide the adjustment slider completely back again, push the volume-setting wheel downward and release the adjustment slider. Remount the cover

NOTICE

The change to factory settings is indicated by the red adjustment slider.

8 *Cleaning and disinfection*

8.1 Autoclaving

The Pipette is completely autoclavable at 121 °C (250 °F), 2 bar and a holding time of at least 15 minutes, in accordance with DIN EN 285.

- a. Eject the pipette tip.
- b. Autoclave the complete pipette without any further disassembling.
- c. Allow the pipette to completely cool and dry.

NOTICE

The effectiveness of autoclaving must be verified by the user. Maximum safety is achieved through vacuum sterilization. We recommend the use of sterilization bags.

If the pipette is autoclaved frequently, the piston and seal should be greased with the supplied silicone grease in order to ensure proper movement. After autoclaving, tighten the connection between the hand grip and the pipette shaft if necessary.

8.2 UV sterilization

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

8.3 PE filter

PE filter for HumaPette 5 ml + 10 ml:

A hydrophobic PE filter is used to prevent liquid from entering the pipette.

Change the filter if it becomes wet or contaminated.

- a. Use a flat object, such as a screwdriver.
- b. Remove the filter carefully, without damaging the tip cone.

Remove the filter before autoclaving!


The device can also be operated without a filter.

9 Maintenance

9.1 Disassembling/cleaning (up to 1,000 μ l)

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks. To do this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms on the pipette tip, refer to *Troubleshooting*, p. 46.

Cleaning

- Detach the pipette shaft (S) from the hand grip by unscrewing it.
- Unscrew the upper part of the ejector unit (A) from the pipette shaft.
- Pull out the shaft (B, C and D) from the lower part of the ejector unit (E).
- Unscrew the piston unit (B).

NOTICE

Do not disassemble the piston unit (B) any further!


- Remove seal with spring (C) (not possible on HumaPette 10 μ l).
- Clean the parts shown with a soap solution or isopropanol, and then rinse with distilled water.
- Dry the parts (max. 120 °C/248 °F).
- Grease piston and seal with a very thin layer of supplied silicone grease.

Reassemble the cooled parts in reverse order. Only hand-tighten the piston unit and the upper part of the ejector unit (A, B).

9.2 Disassembling/cleaning (5 ml/10 ml)

- Check the pipette tip cone for damage.
- Inspect the piston and seal for contamination.
- Check the device for leaks. To do this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms on the pipette tip, refer to *Troubleshooting*, p. 46.

Cleaning

- Remove the entire shaft (S) from the hand grip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
- Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
- Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
- Remove the O-ring-seal from the piston unit and clean it.

NOTICE

Do not disassemble the piston unit (G) any further!


- Clean the piston unit (G) and the bottom part of the shaft (H) with a soap solution or isopropanol, and then rinse with distilled water.
- Dry the parts (max. 120 °C/248 °F) and allow them to cool.
- Carefully lubricate the inside and outside of the O-ring (O) and mount it on the piston.

Reassemble the individual components in reverse order.

10 Troubleshooting

Fault	Possible cause	Corrective action
Tip dripping (device leaking)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Firmly press tip on
The instrument does not aspirate or aspirates too little; the dispensed volume is too low	Seal contaminated	Clean seal
	Seal or cone is damaged	Replace seal or shaft
	Piston is contaminated or damaged	Clean or replace piston
Aspiration is very slow	Shaft is clogged	Clean shaft
	Filter contaminated on 5 ml or 10 ml devices	Change filter
Dispensed volume too large	Pipetting button pressed too far (to the over-stroke point) before aspirating	Ensure proper handling.
Piston sluggish	Piston is contaminated or not greased	Clean piston and apply grease

11 Product markings

Symbol or number	Meaning
	Read the user manual.
	With this mark, we confirm that the product complies with the requirements set out in the EC Directives and has been subjected to the specified testing procedures.
XXZXXXXX	Serial number
	Autoclavable up to the temperature shown

12 Ordering Information

12.1 Order info/accessories

Fixed-volume model

Volume	Description	Order No.
10 μ l	F-10	19505P
20 μ l	F-20	19510P
25 μ l	F-25	19515P
50 μ l	F-50	19520P
100 μ l	F-100	19525P
200 μ l	F-200	19530P
500 μ l	F-500	19540P
1000 μ l	F-1000	19545P

Digital model

Volume	Description	Order No.
0.5 - 10 μ l	D-10	19110P
2 - 20 μ l	D-20	19120P
10 - 100 μ l	D-100	19125P
20 - 200 μ l	D-200	19130P
100 - 1000 μ l	D-1000	19140P
0.5 - 5 ml	D-5000	19150P
1 - 10 ml	D-10000	19160P


Shelf mount for 1 pipette

Description	Order No.
Shelf/rack mount	19945

12.2 Spare parts

12.2.1 Up to 1,000 μ l

The appearance and dimensions of the spare parts correspond to the respective nominal volume. (Fig. HumaPette 20 - 200 μ l)

- A** Ejector (upper part) **B** Piston unit
C Seal with spring **D** Shaft with ejector spring
E Ejector (bottom part)

Fixed-volume model

Volume	A	B	C	D	E
10 μ l	19900-1	19902	–	19914*	19900-4
20 μ l	19900-2	19903	19908	19915	19900-4
25 μ l	19900-2	19907	19913	19915	19900-4
50 μ l	19900-2	19904**		19916	19900-4

Volume	A	B	C	D	E
100 µl	19900-2	19904**		19916	19900-4
200 µl	19900-2	19905**		19917	19900-5
500 µl	19900-3	19906**		19918	19900-6
1000 µl	19900-3	19906**		19918	19900-6

* Seal permanently installed in shaft – not removable!

** Set

Adjustable model

Volume	A	B	C	D	E
0.5 - 10 µl	19900-1	19902	–	19914*	19900-4
2 - 20 µl	19900-2	19903	19908	19915	19900-4
10 - 100 µl	19900-2	19904**		19916	19900-4
20 - 200 µl	19900-2	19905**		19917	19900-5
100 - 1000 µl	19900-3	19906**		19918	19900-6

* 0.1 - 1 µl / 0.1 - 2.5 µl / 0.5 - 10 µl incl. seal

** Set

12.2.2 5 ml/10 ml

The appearance and dimensions of the spare parts correspond to the respective nominal volume. (Fig. replacement parts for HumaPette 5 ml).

Volume	F + F'	G	H	I
0.5 - 5 ml	19934-1	19936	19938	19940
1 - 10 ml	19935-1	19937	19939	19940

12.3 Additional accessories

Description	Order No.
Filter for HumaPette 5 ml, PU 25 pcs.	19941
Filter for HumaPette 10 ml, PU 25 pcs.	19942
Silicone grease for HumaPette up to 1,000 μ l	19880
Silicone grease for HumaPette 5 ml/10 ml	19881

13 Repairs

If any device malfunction cannot be remedied in the lab by simply exchanging replacement parts, the device must be sent in for repair.

Please note that for safety reasons only clean and decontaminated devices can be inspected and repaired!

13.1 Sending for repair

- Clean the device thoroughly and decontaminate.
- Fill out form "Statement on health safety" (pre-printed forms can be requested from dealers or the manufacturer).
- Send the completed form together with the device to the manufacturer or dealer. Include an exact description of the malfunction and the media used.

The device is returned at the risk and expense of the sender.

14 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for

failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

15 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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