


CATALOG NO.....:	150268	
TUBE MATERIAL.....:	POLYSTYRENE	
	INCHES	MM
D1 (OUTSIDE DIA)...:	.69	17.5
D2 (INSIDE DIA)...:	.58	14.7
H1 (HEIGHT).....:	3.746	95.15
NOM. WEIGHT (gms)...:		5.3

CLOSURE MATERIAL...:	LDPE	
	INCHES	MM
A (DIA).....:	.88	22.35
B (HEIGHT).....:	.74	18.8
NOM. WEIGHT (gms)...:		1.3

H2 (TUBE + CLOSURE HEIGHT)...:	3.77	95.8
ASSEMBLED WEIGHT (gms)...:		6.6

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	KOW	APP'D	GRP	PSS 150268	REV. 0711
CK'D	PB		8/4/11		
DESCRIPTION: 15ml TUBE WITH CLOSURE					