



Microplate Washer

Microplate Washer ADX-120

- Programmable soak time, wash cycles, and dispensing volume
- Flexible shaking speed and plate formats
- Auto-clean feature and removeable manifold for easy maintenance
- Continuous aspiration preventing overflow
- Liquid level sensor in waste bottle
- Contains an incubator with two timers
- Specially designed manifold preventing needle from scratching bottom of wells



Specifications

Manifold	8 way manifold autoclavable
Dispensing	Programmable up to 500 µL
Aspirating Method	Continuous through diaphragm pump while dispensing to prevent overflow
Residual Volume	< 2 µL
Waste Bottle	One (capacity 2 L) with audible alarm
Wash Bottle	One (capacity 2 L)
Rinse Bottle	One (capacity 2 L)
Display	16 digit alphanumeric fluorescent
Liquid contact materials	Silicon, stainless steel, derline
Memory	35 open channels
Washing Method	Plate wash, strip wash, bottom wash, overflow wash, rinsing, priming, disinfect
Timer	2 individually programmable timers

ELISA kits



BLOODBORNE

- E0410 HIV 1+2 Ab
- E0420 HIV 1+2 Ag-Ab
- E0511 HCV Ab
- E0610 Syphilis Ab
- E0710 HBsAg



HEPATITIS

- E0100 HAV IgM
- E0105 HEV IgM
- E0510 HCV IgG
- E0511 HCV Ab
- E0810 HBcAb
- E0710 HBsAg
- E0711 HBsAb
- E0720 HBeAg
- E0721 HBeAb
- E0811 HBc IgM



TROPICAL / PARASITIC

- E0170 Chagas IgG
- E0310 Dengue IgM
- E0311 Dengue IgG
- E0312 Dengue Ag
- E0315 CHIK IgM
- E0330 Lepto IgM
- E0320 Zika IgM



FATIGUE / WELLNESS

- E8010 Vitamin D (Total)
- E8030 Vitamin B12
- E1035 TSH
- E5010 Ferritin




HORMONE / NEONATAL

- E1000 hCG
- E1100 FSH
- E1090 LH
- E1080 PRL
- E1010 Total T3
- E1011 Free T3
- E1020 Total T4
- E1021 Free T4
- E1030 TSH
- E1210 Testosterone
- E1220 Progesterone
- E1230 Estradiol (E2)

Distributed by

CTK-MK-FLR-ADX110 Rev 2.1



QUALITY
INNOVATION
SIMPLICITY



Headquarters in California, USA
ctkbiotech.com | info@ctkbiotech.com