



Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

Product Name: Emergency Trolley

Model Number: BIE003A



BIE003A Emergency Trolley:

- High quality ABS and aluminum material
- Mounted on swivel 100mm castors with 2 brakes
- Trolley top with 3 sided rails
- Five drawers with central locking and keys ,Drawers inner with 3x3 divisions
- Side with telescopic working shelf
- Side with bucket, dust buckets and file holder for use
- Top with defibrillator tray and iv pole
- Back with CPR board, electric socket and oxygen tank holder
- Size approx. 750x475x930mm (LxWxH)
- Packed in strong carton

