



Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

Product Name: BED INTENSIVE CARE WITH MATTRESS - ADULT

Model Number: BIH001EA



BIH001EA BED INTENSIVE CARE WITH MATTRESS - ADULT

- Epoxy powder coated steel bed frame, anti rust, strong welding
- Backrest part is radio translucent board for x ray filming
- Mounted on strong and swivel 150mm castors
- Castors are central braking system, with fifth steering castor
- Ergonomic pushing handle on head end
- ABS railing, easily folding down by gas spring system with angle indicator
- Bumpers on all corners of the bed
- Overall dimension: 2230x1050x500-1000mm (LxWxH)
- Heavy loading capacity: 260KG
- Famous denmark Linak motor system
- Controller inserted on bed railing with nurse locking function:
Controller both inside and outside railing for patient and nurse control
- Angle indicator on beds railing part
- Bed comes with night light on both sides of the bed
- Electric CPR and bilateral manual CPR function for emergency rescue
- Bed with rechargeable backup battery
battery can support 80 bed movements
- All beds movements controlled by railing controller:





Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

Backrest lifting: 0-75°

Leg rest lifting: 0-45°

Height adjustment: 500-1000MM

Trendelenburg and reverse Trendelenburg: 0-15°

- Electric One button Cardiac Chair position for patient comfortable
- Bed with fowler and auto contour function
- Electric leg part can be extended 200mm distance
- Bed with 4 hooks telescopic iv pole and bilateral drainage bag hooks
- Class 1, Type B electric shock protection , 220V, 50HZ, 3 pins
- High density foam mattress, good air ventilating for anti bedsores
- Waterproof, fire retardant, anti bacteria, anti puncture, anti static
Anti tearing, anti chemical pu material cover
- Mattress is suitable for bed movements, thickness 200mm
- Packed in strong wood box

BIH001EA BED INTENSIVE CARE WITH MATTRESS - ADULT



X-Ray Translucent Mattress Platform



Tuck away side rail with buffer
Built-in controller and angle indicator



ABS bed ends One Step Removable, meet the first-aid requirements of high-end wards and ICU room



Bumper wheel ,φ74mm, max. protect the bed and patient during transportation



Battery - 18times movement with 100kg bearing, 130times movement without bearing



Full SS IV pole, bearing 15kg





Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

DESIGN

Extra Long and height head foot board, with curve surrounded design, provide comforting feeling for patients



0.12_{mm} | 60° | 50_{kg}



STEELS | PAINTING

Baosteel from Fortune Global 500
|| Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact

WELDING | SAFETY

Panasonic Robotic ensure, 360° full smooth wedding
EN60601 Safety Standard



100_{°C} | 30_{MPa}



PLASTIC

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa

BRAKE

Break spring keep qualified after 10,000times uninteruped break



UL / ROHS / EN

50_{dB} | 20,000_{times}



ELECTRIC SYSTEM

American motor, with UL/ROHS/EN standard. Noise under 50dB, with lifetime 20,000times movement.





Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

