

Speeco Transportation Stretcher

Model:BIPT002H

Hospital
Equipment



Speeco Transportation Stretcher

Solid and accurate performance

Convenient
And
Reliable

Main Features:

- Powerpacker hydraulic pump ensure high quality performance
- One button lay down side rail, gas spring controlled
- PP surface, easy for clean, with mattress could be removable
- Backrest raising with CPR release, dual gas spring controlled
- Central brake sytem at two sides of castors
- Steeringwheel technology, shock-proof
- Hands free trend. & reversed trend.




1 Hands free foot contral pedal at dual sides



2 Central brake pedal at dual ends



3 Directional steering wheel, shock-proof
Foot control pedal at dual sides



4 Dual gas spring control backrest,
with clip for x-ray film



5 Durable button for lay down the
side rail



Optional equip with foldable monitor/
defibrillator tray



Manual crank for knee gatch



Push handle could be folded
down, itegrated IV.pole folded down

Speeco Transportation Stretcher



Standard Features

- PowerPack hydraulic pumps ensure performance
- One button to lay down the side rail
- Hands free hi-lo operation, trend./reversed trend., dual sides control
- Torque release knee gatch handle
- Backrest raising with CPR release
- Anti-slip push handle could be laid down
- Central brake pedal at two ends
- Fifth wheel steering, shock-proof
- Zero sliding gap
- Oxygen cylinder storage, base with tray
- Restraint belt

Standard Accessories

- Stainless steel Telescopic IV pole with four hooks
- Leatherette cover mattress, 8cm thick

Specifications

Overall Length	215cm
Full Width	85cm
Sleeping Deck Length	193cm
Sleeping Deck Width	62cm
Weight Capacity	250kg
Standard Caster	20cm
Fifth wheel	10cm
Max. Height	92cm
Min. Height	60cm

Positioning

Backrest	0°- 90°
Knee Gatch	0°- 40°
Trend./Reverse Trend.	+18°/-18°

BiHealthcare reserves the right to change specifications without notice.

BiHealthcare

MAKE IT SAFER

WARRANTY

- Three years warranty
- 24 hours quick response

ZHANGJIAGANG BRAUN INDUSTRY CO., LTD.



No.24, Zhenxing Road, Zhangjiagang
Economic Development Zone, Jiangsu, China.

Tel: 0086-512-58170517
info@bihealth-care.com

Fax: 0086-512-58170527
web:www.bihealth-care.com

■ As the circumstance of printing, the photos might have minor difference on color, tactile from real products.

■ As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.