

## CMS50D/CMS50D+/CMS50D1/CMS50ED Pulse Oximeter

### Features

- 0.96" dual-color OLED display(blue and yellow).
- Display of SpO<sub>2</sub>, PR, bar graph and pulse waveform.
- Four directions and six display modes.
- Adjustable brightness(only for CMS50D, CMS50D1 and CMS50ED).
- Low battery indication.
- Automatic power off when the finger leaves the device for 5s.
- Display mode can be saved after power off.
- Audible and visible alarm for over-limit(only for CMS50D+).
- PR sound(only for CMS50D+).
- Continuous 24h storage for SpO<sub>2</sub> and PR value, data stored and real-time data can be uploaded to PC(only for CMS50D+).
- Optional enclosure color(only for CMS50D and CMS50D1).
- Power Supply: two "AAA" batteries(only for CMS50D, CMS50D+ and CMS50D1).  
Rechargeable lithium battery, charging time: 2h(only for CMS50ED).
- Continuous working time: CMS50D/CMS50D1: 20h  
CMS50D+: 24h  
CMS50ED: 10h



CMS50D



CMS50D+



CMS50D1



CMS50ED