



CMS6500 Patient Monitor

CMS6500 is a portable monitor with touch screen. To ensure monitoring quality and reduce cooperation risk, it adopts full-isolation (floating), ECG defibrillation-proof protection, ECG anti-high frequency interference and NIBP dual-overpressure protection technique, which is applicable for monitoring adult, pediatric and neonate.

Main Features

- Display mode: 7" color LCD with high-resolution.
- Standard parameters: ECG, SpO₂, NIBP, RESP, TEMP, PR
- Simple and shortcut for operating with touch screen.
- Monitoring parameter, scan speed, volume and output contents can be set
- Optionally. Waveform and parameter color also can be set optionally.
- Storage of 48-hour trend data, and review of 34-second holographic waveform.
- Full-lead ECG displays.
- Anti-high frequency surgical unit, defibrillation-proof(requirement for special leads)
- Digital SpO₂ technology, with strong anti-interference & anti-weak filling capability.
- Networking with Central Nursing station, Bed side and for software updating.
- Built-in rechargeable battery for uninterrupted monitoring.
- Optional functions: Moving bracket, hanging bracket and 12-lead.