

NIBP Series

CONTEC08A/CONTEC08C Electronic Sphygmomanometer

Features

- CONTEC08A : color LCD.
- CONTEC08C : 7-segment LCD.
- Full-automatic NIBP measurement.
- Memory for three-user data, which meets the requirements for whole family.
- Optional units: mmHg and kPa.
- Automatic power-off when there is no operation for several minutes.
- Low battery indication.
- Optional SpO₂ function(SpO₂ probe is necessary).
- Data review interfaces with "Data list"(only for CONTEC08A), "Trend chart"(only for CONTEC08A) and "Big font" ensure parameters are displayed clearly.
- Communicate with PC, to achieve the functions of data review, result analysis, trend chart observation, reports print, etc.
- Dimension: 130mm(L)×110mm(W)×80mm(H)

Only for CONTEC08A

- Measurement mode: adult, pediatric and neonatal.
- Memory: blood pressure results, measurement date and time.

Optional accessories

- Dry battery
- Power adapter
- Ni-MH rechargeable battery
- Ni-MH and Ni-Cd battery charger.
- Neonatal cuff(6-11cm)(only for CONTEC08A)
- Pediatric cuff(10-19cm)(only for CONTEC08A)
- Child cuff(18-26cm)(only for CONTEC08A)
- Big size adult cuff(33-47cm)
- Fingertip SpO₂ probe(adult, child)
- Fingerstall SpO₂ probe(adult)
- Integration SpO₂ probe
- Disposable SpO₂ probe(adult, neonatal, child and pediatric)

