

CONTEC08D/CONTEC08E Electronic Sphygmomanometer

Features

- CONTEC08D : 7-segment LCD.
- CONTEC08E : color LCD.
- Be applicable for adult.
- Start measuring manually, record each measurement data, and store up to 99 groups of data.
- Low battery and error information indication.
- Optional SpO₂ function: switch to main interface automatically after inserting the SpO₂ probe.
- Automatic power-off when there is no operation for several minutes.
- Optional units: mmHg and kPa.
- Memory for blood pressure results, measurement date and time, BP list and trend chart with small font(only for CONTEC08E).
- Dimension: 132mm(L)×130mm(W)×50mm(H)

Optional accessories

- Big size adult cuff(33-47cm)
- Fingertip SpO₂ probe(adult, child)
- Fingerstall SpO₂ probe(adult)
- Integration SpO₂ probe
- Disposable SpO₂ probe(adult, neonatal, child and pediatric)



CONTEC08D



CONTEC08E