

BESDATA

BESDATA® BD-DF Reusable Video Laryngoscope



1

The ergonomic design, easy to intubation, easy to hold.

2

6 size blades could support with Routine and Difficult Airways Intubation.

3

Compatible with BD-M3 3inch portable monitor and DB-M7 7inch separate big monitor.

4

More Available Space with 12mm BESDATA blades height

5

Anti-fog blade and camera.



Model Number: BD-DF

Order Number	BD-M3	BD-DF-R5	BD-DF-R4	BD-DF-R3	BD-DF-R2	BD-DF-R1
Description	3inch LCD Monitor	Miller 0	MAC1	MAC2	MAC3	MAC4
Application	LCD Monitor	Neonate	Pediatric	Small Adult	Adult	Large Adult
Blade Length	/	8.2	10.45	12.11	14.25	16.08

Technical Data

Parts	Description	Technical Index
Monitor	Display	3 inch LCD
	Resolution Ratio	640*480(RGB)
	Viewing angle	>60°
	Power	<2W
	angle	0°~130°
	Display up and down rotation angle	0-270°
Camera	Illuminance	>150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	/video
	Memory type	Micro SD card 8GB
	Video/Photo Play back	Yes
Batterty	Battery Type	Rechargeable lithium battery
	Capacity	1350mAh
	Cycle life of battery	>300 times
	Battery working time	>140min
	Charging time	<150min