

WARRANTY

- Three years warranty
- 24 hours quick response

Electric bed 5 function / 3 function

Model:BIH007EB-5/BIH007EB-3

5 KEY ADVANTAGES

- DESIGN & QUALITY
- COST-EFFECTIVE PRICE
- HIGHLY SUPPORT ON CUSTOMIZED PRODUCTS
- 3-YEAR SPARE PARTS WARRANTY
- WHOLE-LIFE TIME FRAME & BASE



Standard Features

- 3 motors for 3-function/4 motors for 5-function
- Hand controller with backlight
- Detachable head and foot board with bumper
- Aluminum alloy side rail, one hand control to lay down
- Four section ventilate platform, epoxy coated steel
- Epoxy coated steel frame and base
- Luxury covered castors with individual brake
- With drainage hook and mattress holder
- With IV. pole sockets at four corners

Optional Feature


- Manual CPR quick release handle at both sides
- Back-up battery

Standard Accessories


- Stainless steel Telescopic IV pole with four hooks

Specifications:

Overall length: 215cm
 Overall width: 96cm
 Caster diameter: 12.5cm
 Bed platform length: 194cm
 Bed platform width: 90cm
 Side rail: 148(L)*41.5(H)cm
 Siderail height over platform: 32cm
 Height range:
 Max.Height: 75cm
 Min.Height: 45cm
 Positions:
 Backrest raising: 0-75°
 Knee raising: 0-45°
 Trendelenburg: 0-12°(5-function available)
 Reverse trendelenburg: 0-12°(5-function available)



5-function panel



3-function panel



Handle with backlight



Push red button to lay down down the side rail

ZHANGJIAGANG BRAUN INDUSTRY CO., LTD.



Buliding C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang City, Suzhou City, Jiangsu Province, China

Tel: 0086-512-58170517
info@bihealth-care.com

Fax: 0086-512-58170527
web:www.bihealth-care.com

■ As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
 ■ As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.