



# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

**Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..**

**Product Name: Electric Bed Complete**

**Model Number: BIH008EA**



## **BIH008EA Electric Bed Complete**

- Epoxy powder coated steel bed frame, anti rust, strong welding
- Four sections perforated steel bed platform with mattress retainers
- Mounted on strong and swivel 125mm castors
- Castors are central braking system, with steering
- ABS railing, easily folding down by gas spring system with angle indicator
- Bumpers on all corners of the bed
- Overall dimension: 2180x1050x450-750mm ( LxWxH )
- Heavy loading capacity: 260KG
- Strong and low noise Germany Dewert motors
- Controller inserted on bed railing with nurse locking function:  
Controller both inside and outside railing for patient and nurse control
- Hand hold wired remote controller for use
- CPR function for emergency rescue, detachable headboard for CPR
- Bed with rechargeable backup battery  
battery can support 80 bed movements
- All beds movements controlled by railing controller&hand hold controller:  
Backrest lifting: 0-75°  
Legrest lifting: 0-35°  
Height adjustment: 450-750MM





# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)


Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

**Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..**

Trendelenburg and reverse trendelenburg: 0-15°

- Bed with fowler and auto contour function
- Bed with 4 hooks telescopic iv pole and bilateral drainage bag hooks
- Class 1, Type B electric shock protection , 220V, 50HZ, 3 pins
- Waterproof, fire retardant, anti bacteria, anti puncture, anti static  
Anti tearing, anti chemical pu material cover
- High density foam mattress, thickness 120mm
- Packed in strong wood box

## BIH008EA Electric Bed Complete





# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

**Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..**

Outside panel




Inside panel



Embedded controls for nurse and patient

foot end



both sides



Soft contoured mattress retainer handles



Has remote control aside from the buttons at the side of the bed





# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

**Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..**

