


Electric 5-function/3-function bed

Model:BIH007EC

Hospital
Equipment



5-function handle



3-function handle

Specifications:

Overall length: 212cm
Overall width: 105cm
Caster diameter: 12.5cm

Height range:
Max.Height: 70cm
Min.Height: 40cm

Positions:
Backrest raising: 0-75°
Knee raising: 0-35°

Trendelenburg: 0-12° (5-function available)
Reverse trendelenburg: 0-12° (5-function available)

Standard Features

- Actuator control for movements, with hand controller
- Detachable head and foot board with bumper
- Split side rail, one handed control to lay down
- Four section ventilate platform, epoxy coated steel
- Epoxy coated steel frame and base
- Central brake double sides castors
- With drainage hook
- With Iv pole sockets at four corners

Standard Accessories

- Stainless steel Telescopic IV pole with four hooks

BiHealthcare

MAKE IT SAFER

WARRANTY

- Three years warranty
- 24 hours quick response

ZHANGJIAGANG BRAUN INDUSTRY CO., LTD.

Building C7, No.199,Beihuan Road,
Zhangjiagang,Jiangsu 215600 P.R.China.

Tel: 0086-512-58170517
info@bihealth-care.com

Fax: 0086-512-58170527
web:www.bihealth-care.com

■ As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
■ As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.