

Stainless Steel Stretcher

Model:BIPT006S



Standard Features

- Whole stretcher body is made by stainless steel
- Mechanical structure for the hi-low adjustment
- Stretcher surface could be separated for emergency
- One button to lay down side rail
- With utility tray and oxygen holder
- Solid castors with brakes
- Comfortable mattress fixed with surface
- IV stand could storage under the surface

Standard Accessories

- Stainless steel Telescopic IV pole with two hooks
- Leatherette cover mattress, thich 6cm

Specifications:

Overall length: 212cm
Overall width: 80cm

Height range:
Max.Height: 85cm
Min.Height: 60cm

 BiHealthcare

WARRANTY

- Three years warranty
- 24 hours quick response

ZHANGJIAGANG BRAUN INDUSTRY CO., LTD.

No.24, Zhenxing Road, Zhangjiagang
Economic Development Zone, Jiangsu, China.

Tel: 0086-512-58170517
info@bihealth-care.com



Fax: 0086-512-58170527
web:www.bihealth-care.com

■ As the circumstance of printing, the photos might have minor difference on color, tactile from real products.