



# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No 199, Hongwu Avenue, Tangqiao Town, Zhangjiagang City, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [info@bihealth-care.com](mailto:info@bihealth-care.com) [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

**Product Name: IV Stand**

**Model Number: BIHR001B**



## **BIHR001B IV Stand:**

- High quality chromed steel frame
- Top with 4 hooks
- Hooks is heavy loading and easy folding
- Base mounted on swivel 50mm castors with brakes
- Stable and strong loading stainless steel base
- Height adjustable by manual knob: 1200-2000mm
- Packed in strong carton

