



Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No 199, Hongwu Avenue, Tangqiao Town, Zhangjiagang City, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: info@bihealth-care.com daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..



BIH010MG Hospital Bed:

- Epoxy powder coated steel frame, anti rust, easy cleaning
- Detachable ABS head and foot boards with bumpers
- Split ABS railing, with angle indicator
- Mounted on swivel 125mm castors with central braking
- Four sections perforated steel bed top with mattress retainers
- Bed movements controlled by foot end retractable cranks:

Backrest lifting: 0-80, Legrest lifting: 0-40

Height adjustment: 450-750mm

- Overall size approx.: 2180x1050x450-750mm(LxWxH)
- Bed with iv pole 4 hooks and drainage bag hooks
- Packed in strong wood box

Product Name: Hospital Bed

Model Number: BIH007EC

