



Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

Product Name: Operating Table

Model Number: BIOT8000



BIOT8000 Operating Table:

- Longitudinal slide up to 350mm for free access to C arm and X ray
- High quality 304 grade stainless steel base , column
- Weight load capacity of 280kgs with highest safety and stability.
- Memory pad with antistatic, waterproof design.
- Standard override panel control and remote controller
- Split leg section for different surgery
- Built-In Kidney Bridge, More Convenient And Nice Appearance.
- Pedal Brake, Easy To Move And Lock.
- Extra-Long Horizontal Sliding, Available For C-Arm Operation.
- All Function Controlled By Hand remote and column Controller
- Mounted on swivel 150mm castors with central manual brakes
- External Size (LxWxH) 2060*500*(710-1010)mm
- Reverse Trendelenburg 0-25° ,Trendelenburg 0-25°
- Lateral Tilt Left/right 0-15°
- Head Plate up/down 0-40° /0-90°
- Back Plate up/down 0-80° /0-20°
- Leg Plate up/down 0-15° /0-90°
- Kidney Bridge:110mm, backup battery in case of power fail
- Standard anesthesia frame, shoulder support, arm support, body support, leg support, iv stand, abdominal safety strap, collection bowl
- Packed in strong wood box

