

Ovapit Handling Pipettes

OVApit is a Manipulation pipette used to handle Huge cells in ART procedures.
Supplied in a peel-open package.

- **Made of flexible Polypropelyne**
Designe to be flexible, not to scratch the dish surface
- **Intensive tests criteria for safer gametes and embryo handling**
Each pipettes are inspected under extreme magnification, sterilized by gamma irradiation andn lot tested for one-cell MEA survivability and LAL toxicity.
- **Each pipette contains a syringe as a hanfling device and the hadling tip to be connected to it for easy of use**

