

Table Top Centrifuge

PLC-01~ 05

GMDN Code 36465

- **The Best-Selling & Most-Valued Product**
- **Cheap Price & Excellent Quality, Guaranteed**
- **2 types of motors for choice: 67W, 90W**
- **For general purpose and also PRP**

Intended use :

A device that is a general-purpose laboratory centrifuge used to separate the components of suspensions by the application of centrifugal force. The device is intended to centrifuge patient samples, e.g. body fluids, either alone or after addition of reagents or other additives before measuring analytes. Ideal for use in hospital, clinics, blood bank, laboratory center, school, veterinary, biochemistry, industry, etc.



Specification

Code	PLC series
Dimension	26 x 26 x 26.5cm approx.
Rotor	Fixed angle
Timer (Mechanical)	60 minutes
Speed	1000 ~ maximum, adjustable
Safety shutout system	Option (shutout + Locking lid)
Standard accessory	Tube adaptors
Standards	ISO9001, ISO13485, GMP,, CE, USFDA
Net weight	8.0 kg or 17.64 lbs.

Option :

- Safety Locking Lid
- AC / DC use in car
- Motor 67W or 90W



Tube Adaptor: (Dimension)

- 17.5mm (I.D.) x 88mm (I.L): For 8.5ml ~ 15ml tubes (21.5mm (O.D.) x 94mm (Outer Length))
- 14mm (I.D.) x 66mm (H): For 3ml ~ 7ml tubes
- Insert 14mm adaptor into 17mm adaptor for more use

Ordering Information

Model	Capacity	Rotor	Max rpm	Max RCF
PLC-01	4 x 15ml	A-0415	4500	2422xg
PLC-02	6 x 15ml	A-0615	4500	2422
PLC-03	8 x 15ml	A-0815	4500	2422
PLC-04	10 x 15ml	A-1015	4000	2057
PLC-05	12 x 15ml	A-1215	4000	2057

Useful Adaptor



For PRP








- ← 3 ~ 7ml Vacutainer
- ← 14 Ø x 66mm adaptor
- ← 17.5 Ø x 88mm adaptor (Outer: 21.5 x 94mm)

Technical Data

- Speed regulated from 1000rpm to the maximum within 15 seconds.
- Choice of variable speed regulation or fixed speed at 3000 or anywhere as requested.
- Anti-vibration rubber buffer balances rotor when spinning, and also protect equipment during transportation.
- Timer activated centrifuge, the LED pilot lamp indicates the power is switched on.
- The heavy casting iron base and 3 rubber suction feet prevent the unit from creeping and moving.
- This unit complies with CE electromagnetic compatibility requirement of IEC/EN61010-1-1:2010 and council medical directive 2004/108/EC.

Tube number available for each rotor?

Volume		10~15ml		10~12ml		8.5~10ml		4~8ml		4~8ml
Dimension		17x120		17x100		16x100		13x66		13x100
PLC-01		4 tubes		4 tubes		4 tubes		4 tubes		4 tubes
PLC-02		2 tubes		6 tubes		6 tubes		6 tubes		6 tubes
PLC-03		4 tubes		8 tubes		8 tubes		8 tubes		8 tubes
PLC-04		4 tubes		10 tubes		10 tubes		10 tubes		10 tubes
PLC-05		4 tubes		12 tubes		12 tubes		12 tubes		12 tubes
Adaptor		17.5 x 88mm		17.5 x 88mm		17.5 x 88mm		14 x 66mm		17.5 x 88mm