



Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need.

Product Name: Reclinder

Model Number: BIHA018A



BIHA018A Reclinder:

- High quality epoxy powder coated steel frame
- Ergonomic armrest for use
- Mounted on swivel castors with brakes
- Electric adjusted back, leg and Trendelenburg position
- Back with pillow , leg section with foot support
- Seat is approx 60x60cm, height is 60cm
- Back is 60cm width , approx 70cm height
- With CPR quick release , Trendelenburg function
- Recliner with side tray
- Strong and stable to use
- Packed in strong carton

