



GEMMY INDUSTRIAL CORP.

High Pressure Autoclave Sterilizer



SA-232



TC



SA-300H



SA-232X



HL-340, 341



SA-252F



HL-342, 343



SA-300VF



SA-300MA



SA-260FA, 260MA

Specifications

Model	SA-232	SA-232X	SA-252F	SA-300H	SA-300VF	SA-300VM	SA-230MA SA-260MA	SA-230FA SA-260FA
Capacity	16 L.	16 L.	24 L.	40 L.	50 L.	40 L. Class B	16 L. Class B 24 L. Class B	16 L. 24 L.
Temperature (°C)	122°C only	115 – 134	121 – 132	121 – 132	121 – 132	121/135°C	121 · 135 121 · 135	121 · 135 121 · 135
Drying	Manual	Manual	Automatic	Automatic	Automatic	Automatic	Automatic	Automatic

Model	HL-340 / 341	HL-342 / 343	HL-342P	HL-343P	TC series (All stainless steel)
Capacity	49 L / 72 L. Stainless steel	49 L / 72 L.	49 L	72 L	530; 424; 339; 283; 196 143; 113; 83; 77; 56 L.
Temp (°C)	120 – 132°C	120 – 132°C	120 – 132°C	120 – 132°C	121 – 134°C
Auto Drying	Option	Option	Automatic	Automatic	Option

- Ideal for sterilization of wrapped, unwrapped, liquid, surgical tools, metal, labwares etc. in hospital and clinics.
- TC series functions automatic drying, safety lid lock by pressure, automatic water filling, pre-vacuum & Dry vacuum etc. which are optional for extra costs.
- SA-300VM, SA-230MA, SA-260MA are Class B autoclaves.
- Quality guaranteed and services guaranteed.
- Contact us for more details.

