

CONTEC SP70B handheld portable medical BT spirometer

Product Description

CONTEC™



SPIROMETER

is a hand-held equipment for checking lung conditions, it adopts the infrared mode for measuring relative items, and is applicable for hospital, clinic, family for routine test.

Features



- Parameters including FVC 、 PEF 、 FEV1 、 FEV1%、 FEF25、 FEF50、 FEF75 and FEF2575 can be measured.
- Personal information setup, including height, age, gender, etc.
- Health status indication.
- Bluetooth transmission.
- Data memory and case review function.
- Low battery indication.
- Real-time clock display and setup.
- Automatic power off function.
- Mass storage device, stores up to 9999 cases.

Details show



Performance



- Volume range: 0%10 L
- Flow rate range: 0 L/s%16 L/s
- Volume accuracy: $\pm 3\%$ or 0.05 L (whichever is greater)
- Flow rate accuracy: $\pm 5\%$ or 0.2 L/s (whichever is greater)
- Safety classification: internally powered equipment, type BF applied part
- Working current: 80 mA
- Power supply: two AAA batteries