



# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

**Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..**

**Product Name: Baby Cot**

**Model Number: BIHB003A**



## **BIHB003A Baby Cot:**

Mobile baby cot with removable bassinet.  
Bassinet is transparent for easy watching  
Bassinet with holes for good ventilation  
Mounted on 4 swivel anti-static castors 75mm diameter  
Swivel and non-mark castors with 2 brakes  
Padded mattress, detachable for easy cleaning.  
Mattress cover removable via side zipper  
Basinet sets and removes smoothly from cart frame.  
High resistance to corrosion  
Epoxy oven baked powder coated tubular steel frame  
High-density polyurethane foam with density approx. 30 kg/m<sup>3</sup>.  
Cover: flexible highly tear resistant, anti-static, flame retardant,  
disinfectant- and liquid proof, washable.  
Bassinet Size: 800x450x250mm ( LxWxH )  
Overall cot size approx: 820x460x900mm ( LxWxH )  
Knockdown design, easy installation  
Stable connection bar on lower part to ensure stability  
Packed in strong carton

