



# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

**Product Name: Blood Donor Chair**

**Model Number: BIHA410A**



## **BIHA410A Blood Donor Chair:**

- ● Multi-position adjustment, using high-performance imported silent push rod motor, easily adjust the backrest, leg rest, footrest, one-button CPR, one-button reset.
- The buttons of hand control is simple, intuitive and easy to operate.
- International brand silent 24V dc push rod motor, work stably, safe
- Long - term continuous working system.
- 10 years of life design, reduce the cost of maintenance and replacement.
- Ultra-low energy consumption, using less than 0.12 degrees per day.
- The chair cushion is made of high density sponge, the elasticity is moderate, sitting and lying for a long time will not cause postural strain; PVC leather is durable, luxurious and comfortable.
- Four motors to adjust backrest,legrest,seat height,footrest  
Backrest adjustment :-10° ~75°  
Seat height adjustment:580-800mm  
Leg rest adjustment :-70° ~10°  
Trendelenburg: -10°
- Total Length: 2100mm Seat Width:600mm(930mm include armrests)
- with:IV POLE ,Paper roll holder
- Mounted on swivel castors central braking pedal
- Packed in strong wood box

