

## BESDATA® BD-M7DF Reusable Video Laryngoscope with 7 Inch Monitor



1

Reusable , Durable Material

2

Anti-Fog Design

3

Separate 7inch LCD Monitor,  
to keep doctor away from  
patients secretions.

4

More Available Space with  
12mm BESDATA blades height

5

Blade and Handle assists with  
either Routine and Difficult  
Airways

6

Blade work both with BD-M3  
portable monitor and BD-M7  
7inch monitor



Model Number: BD-M7DF

Order number	BD-M7	BD-DF-R5	BD-DF-R4	BD-DF-R3	BD-DF-R2	BD-DF-R1	BD-DF-R6
Discription	7inch LCD Monitor	Mil0	MAC1	MAC2	MAC3	MAC4	D Blade
Application	LCD Monitor	Neonate	Pediatric	Small Adult	Adult	Big Adult	Difficult Airway
Blade Length	/	8.2cm	10.5cm	12.1cm	14.2cm	16.1cm	13.2cm

## Technical Data

Parts	Description	Technical Index
Monitor	Display	3 inch LCD
	Resolution Ratio	640*480(RGB)
	Viewing angle	>60°
	Power	<2W
	Display front and back rotation angle	0°~130°
	Monitor Thickness	12mm
	Weight	<90g
Camera	Illuminance	>150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	USB output, Easy to storage and establish image /video
	Memory type	Micro SD card 8GB
	Video/Photo Play back	Yes
Battery	Battery Type	Rechargeable lithium battery
	Capacity	730mAh
	Cycle life of battery	>300 times
	Battery working time	>75min
	Charging time	<90min
Blades	Weight	< 80g