



# Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No 199, Hongwu Avenue, Tangqiao Town, Zhangjiagang City, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: [info@bihealth-care.com](mailto:info@bihealth-care.com) [daniel@bihealth-care.com](mailto:daniel@bihealth-care.com)

Skype: dan.bihealthcare Whatsapp ( Wechat, Viber ): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need..

**Product Name: Examination Table**

**Model Number: BIDE003A**



## **BIDE003A Examination Table:**

- High quality epoxy powder coated steel frame.
- Strong welding, anti rust, easy cleaning
- Waterproof, fire retardant pu cover high density foam table top
- Two sections bed platform
- Back section lifted by manual rack
- Paper roll holder for use
- Rubber protective caps on four legs
- Heavy loading capacity
- Approx size: 1900x600x700mm(LxWxH)
- Packed in strong carton

