



Zhangjiagang Braun Industry Co.,Ltd.

Building C7, No.199, Hongwu Avenue, Tangqiao Town, Zhangjiagang, Jiangsu Province, China.

Tel: 86 512 58170517 Fax: 86 512 58170527

Email: daniel@bihealth-care.com

Skype: dan.bihealthcare Whatsapp (Wechat, Viber): 86 18626299207

Focus on Medical Supply, Provide Best Medical Solution. We Care What you Need.

Product Name: Examination Couch

Model Number: BIDE003A



BIDE003A Examination Couch:

- High quality epoxy powder coated steel frame
- Strong welding with no sharp edges
- Two sections couch top
- Top is waterproof, fire retardant pu cover high density foam
- Backrest is lifting by manual rack
- Approx size: 2000x600x680mm(LxWxH)
- Paper roll holder for use
- Heavy loading capacity
- Durable and strong to use
- Packed in strong carton

