





Safepit denudation Pipettes

	Size of ID (um)	Catalogue Number	Intended Use
	200	Size of ID (um)	Balstocyst & And Embryo Handling
	300	Size of ID (um)	Balstocyst & And Embryo Handling
	170	Size of ID (um)	Handling Embryos
	145	Size of ID (um)	Oocyte Stripping For ICSI

Safepit is a Manipulation pipette used to handle cells and remove cumulus cells in ART procedures. Supplied in a peel-open package.

- **Made of flexible Polycarbonate**
Designe to be flexible, not to scratch the dish surface
- **Intensive tests criteria for safer gametes and embryo handling**
Each pipettes are inspected under extreme magnification, sterilized by gamma irradiation and lot tested for one-cell MEA survivability and LAL toxicity.
- **Universal size that could be used with any handling device.**